

## **n-Decane - Flash point**

Art. ID ERM-FC033  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2247 III E1

### **Description**

The certified value was determined by the Abel closed cup method described in the Institute of Petroleum Standard IP170/95 and also published as British Standard BS2000:Part 170: 1995. The certified value is corrected to standard barometric pressure at 0°C.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Non-equilibrium flash point		32,5 ± 0,5	°C		