

Gas oil - Cold filter plugging point (CFPP) and cloud point (CP)

Art. ID	ERM-FC395K
Unit	set(2 x 27 mL)
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Description

Intended for verifying manual and automatic cold filter plugging point (CFPP) apparatus in relation to the methods EN 116 / ASTM D6371 and cloud point (CP) in relation to the methods ISO 30115 and ASTM D2500.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cold filter plugging point (CFPP)		-7,9	°C	As defined by EN 116 or ASTM D 6371 using the automatic or manual procedure	
	Cloud point (CP)		-7,2	°C	As defined by EN 230 15, ISO 30 15 or ASTM D2500 using the automatic or manual procedure	