

Acrylonitrile-d3, 98 atom % D, >=99% (CP), contains hydroquinone as stabilizer (ALDRICH)

Art. ID SAF-485314-1G
Unit 1 g
Deliverydetails ADR Dangerous Good (DG) / AIR Cargo Air Only (CAO) UN1093 (6.1) I E0

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
| | Acrylonitrile-d3 | [53807-26-4] | | | | |