

**Sodium L-lactate-13C3 solution,45-55 % (w/w) in H2O, >=99 atom % 13C, >=98% (CP), >=98% (Chiral purity, HPLC) (ALDRICH)**

Art. ID                      SAF-485926-500MG  
Unit                         500 mg  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sodium L-lactate-13C3	[201595-71-3]	45 - 55	%		