

Tin Base, CRM, Solder, Disc 40 mm

Art. ID IMN-L-3
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Arsenic (As) | [7440-38-2] | 0,092 | % | | |
| | Bismuth (Bi) | [7440-69-9] | 0,22 | % | | |
| | Cadmium (Cd) | [7440-43-9] | 0,0065 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,034 | % | | |
| | Iron (Fe) | [7439-89-6] | ~0,023 | % | | |
| | Antimony (Sb) | [7440-36-0] | 0,14 | % | | |
| | Tin (Sn) | [7440-31-5] | 60,18 | % | | |
| | Zinc (Zn) | [7440-66-6] | 0,0064 | % | | |