

**Sodium L-lactate-3-13C solution,45-55 % (w/w) in H2O, >=99 atom % 13C, >=98% (CP), >=98% (Chiral Purity, HPLC) (ALDRICH)**

Art. ID                      SAF-490040-250MG  
Unit                         250 mg  
Deliverydetails            No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter      | CAS number    | Concentration/Value | Unit | Method | Source |
|------------------|------------------------|---------------|---------------------|------|--------|--------|
|                  | Sodium L-lactate-3-13C | [201595-70-2] | 45 - 55             | %    |        |        |