

Sodium L-lactate-3-13C solution,45-55 % (w/w) in H2O, >=99 atom % 13C, >=98% (CP), >=98% (Chiral Purity, HPLC) (ALDRICH)

Art. ID SAF-490040-250MG
Unit 250 mg
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|------------------------|---------------|---------------------|------|--------|--------|
| | Sodium L-lactate-3-13C | [201595-70-2] | 45 - 55 | % | | |