

1-Iodopropane-d7, >=98 atom % D, >=99% (CP), contains copper as stabilizer (ALDRICH)

Art. ID SAF-614734-5G
Unit 5 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2392 III

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
| | 1-Iodopropane-d7 | [59012-23-6] | | | | |