

**Sodium L-lactate-3,3,3-d3 solution,45-55 % (w/w) in H2O, >=98 atom % D, >=98% (CP), >=98% (Chiral Purity, HPLC) (ALDRICH)**

Art. ID                      SAF-616702-1G  
Unit                         1 g  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sodium L-lactate-3,3,3-d3		45 - 55	%		