

Silver Standard for ICP-MS, TraceCERT®, 1 mg/L Ag in nitric acid (SIGMA-ALDRICH)

Art. ID SAF-69389-100ML
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

limited shelf life, expiry date on the label

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	1	mg/L		