

**Acrylamide-d3 Standard solution, ~500 mg/L in acetonitrile, analytical standard (SIGMA-ALDRICH)**

Art. ID                      SAF-72834-5ML-F  
Unit                         5 mL  
Deliverydetails            ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1648 II E2

| Text/Information | Analyte/Parameter | CAS number    | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|---------------|---------------------|------|--------|--------|
|                  | Acrylamide-d3     | [122775-19-3] | ~500                | mg/L |        |        |