

**Acrylamide-d3 Standard solution, ~500 mg/L in acetonitrile, analytical standard (SIGMA-ALDRICH)**

Art. ID                      SAF-72834-5ML-F  
Unit                         5 mL  
Deliverydetails            ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1648 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acrylamide-d3	[122775-19-3]	~500	mg/L		