

**Tetrahydrocortisone-2,2,3,4,4-d5 solution,100 µg/mL in methanol, >=98 atom % D, >=98% (CP) (ALDRICH)**

Art. ID                      SAF-809276-1ML  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tetrahydrocortisone-2,2 ,3,4,4-d5		100	µg/mL		