

Lead in gasoline - 4 g/gal

Art. ID ASI-PBG7A-06-100ML

Unit 100 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1203 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|------------------------|--------|
| | Lead (Pb) | [7439-92-1] | 4 | g/gal | ASTM D5059 , PART A | |