

Lead in isooctane - 0.3 g/gal

Art. ID ASI-PBIO7C-07-100ML

Unit 100 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1262 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|------------------------|--------|
| | Lead (Pb) | [7439-92-1] | 0,3 | g/gal | ASTM D5059 , PART A | |