

Nickel Reagent and 0.05% Manganese internal standard in xylene-mineral oil - Ni 0.005%

Art. ID ASI-NiISO7-05-100ML
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1307 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|------------|--------|
| | Nickel (Ni) | [7440-02-0] | 0,005 | % | ISO/CD 145 | 96 |