

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Aluminum 3%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-AI-3.0-25G  
Unit 25 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	3	%		