

Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Chromium 3.5%, Sulfur-Free < 1 ppm

Art. ID ASI-LOMS-Cr-3.50-25G
Unit 25 g
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Chromium (Cr) | [7440-47-3] | 3,5 | % | | |