

Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Lithium 1.5%, Sulfur-Free < 1 ppm

Art. ID ASI-LOMS-Li-1.5-25G
Unit 25 g
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lithium (Li)	[7439-93-2]	1,5	%		