

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) -  
Manganese 6%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-Mn-6.0-25G  
Unit 25 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Manganese (Mn)	[7439-96-5]	6	%		