

Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Strontium 10%, Sulfur-Free < 1 ppm

Art. ID ASI-LOMS-Sr-10.0-25G
Unit 25 g
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Strontium (Sr)	[7440-24-6]	10	%		