

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Zirconium 5%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-Zr-5.0-25G  
Unit 25 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Zirconium (Zr)	[7440-67-7]	5	%		