

Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Barium 12.5%, Sulfur-Free < 1 ppm

Art. ID ASI-LOMS-Ba-12.5-50G
Unit 50 g
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Barium (Ba) | [7440-39-3] | 12,5 | % | | |