

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Cerium 5%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-Ce-5.0-50G  
Unit 50 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Cerium (Ce)       | [7440-45-1] | 5                   | %    |        |        |