

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Iron 4%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-Fe-4.0-50G  
Unit 50 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Iron (Fe)	[7439-89-6]	4	%		