

Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Magnesium 3%, Sulfur-Free < 1 ppm

Art. ID ASI-LOMS-Mg-3.0-50G
Unit 50 g
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Magnesium (Mg) | [7439-95-4] | 3 | % | | |