

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Silicon 7.5%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-Si-7.5-50G  
Unit 50 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	7,5	%		