

Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Vanadium 4%, Sulfur-Free < 1 ppm

Art. ID ASI-LOMS-V-4.0-50G
Unit 50 g
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Vanadium (V) | [7440-62-2] | 4 | % | | |