

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Zinc 6%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-Zn-6.0-50G  
Unit 50 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Zinc (Zn)	[7440-66-6]	6	%		