

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Cerium 5%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-Ce-5.0-100G  
Unit 100 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cerium (Ce)	[7440-45-1]	5	%		