

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Chromium 3.5%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-Cr-3.50-100G  
Unit 100 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chromium (Cr)	[7440-47-3]	3,5	%		