

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Copper 6%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-Cu-6.0-100G  
Unit 100 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	6	%		