

Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Iron 4%, Sulfur-Free < 1 ppm

Art. ID ASI-LOMS-Fe-4.0-100G
Unit 100 g
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Iron (Fe) | [7439-89-6] | 4 | % | | |