

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Lead 20%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-Pb-20.0-100G  
Unit 100 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lead (Pb)	[7439-92-1]	20	%		