

Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Magnesium 3%, Sulfur-Free < 1 ppm

Art. ID ASI-LOMS-Mg-3.0-100G
Unit 100 g
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Magnesium (Mg)	[7439-95-4]	3	%		