

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Potassium 7.5%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-K-7.5-100G  
Unit 100 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Potassium (K)     | [7440-09-7] | 7,5                 | %    |        |        |