

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Selenium 3.5%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-Se-3.5-100G  
Unit 100 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Selenium (Se)	[7782-49-2]	3,5	%		