

**Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Sodium 2.5%, Sulfur-Free < 1 ppm**

Art. ID ASI-LOMS-Na-2.5-100G  
Unit 100 g  
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sodium (Na)	[7440-23-5]	2,5	%		