

Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Vanadium 4%, Sulfur-Free < 1 ppm

Art. ID ASI-LOMS-V-4.0-100G
Unit 100 g
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Vanadium (V)	[7440-62-2]	4	%		