

Concentrated Liquid Organo-Metallic Solutions (Stabilized and Ready for Use.) - Yttrium 2.5%, Sulfur-Free < 1 ppm

Art. ID ASI-LOMS-Y-2.5-100G
Unit 100 g
Deliverydetails /not restricted

Description

These concentrated solutions are stabilized and can be used to formulate sets of standards for the analysis of additive elements in lubricating oils

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Yttrium (Y)	[7440-65-5]	2,5	%		