

Sulfur in Diesel Fuel Calibration Standards, Conc. Range :0 - 5.0 wt%, 10 standards per set, 100 mL each (Standard for XRF)

Art. ID ASI-SDF10-Set(10x100ML)
Unit 10 x 100 mL
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Standard 1	Sulfur (S)	[7704-34-9]	0	%	ASTM D2622 & D4294	
Standard 2	Sulfur (S)	[7704-34-9]	0,1	%	ASTM D2622 & D4294	
Standard 3	Sulfur (S)	[7704-34-9]	0,25	%	ASTM D2622 & D4294	
Standard 4	Sulfur (S)	[7704-34-9]	0,5	%	ASTM D2622 & D4294	
Standard 5	Sulfur (S)	[7704-34-9]	0,75	%	ASTM D2622 & D4294	
Standard 6	Sulfur (S)	[7704-34-9]	1	%	ASTM D2622 & D4294	
Standard 7	Sulfur (S)	[7704-34-9]	2	%	ASTM D2622 & D4294	
Standard 8	Sulfur (S)	[7704-34-9]	3	%	ASTM D2622 & D4294	
Standard 9	Sulfur (S)	[7704-34-9]	4	%	ASTM D2622 & D4294	
Standard 10	Sulfur (S)	[7704-34-9]	5	%	ASTM D2622 & D4294	