

Sulfur (S) (as Di-N-butyl sulfide) in residual oil. Conc. Range: 0.25 - 5.00 wt%, 10 standards per set, 100 mL each (Standard for XRF)

Art. ID ASI-SRO10-Set(10x100ML)
Unit 10 x 100 mL
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-SRO10-01-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,25	%	ASTM D2622 & D4294	
ASI-SRO10-02-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,5	%	ASTM D2622 & D4294	
ASI-SRO10-03-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	1	%	ASTM D2622 & D4294	
ASI-SRO10-04-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	1,5	%	ASTM D2622 & D4294	
ASI-SRO10-05-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	2	%	ASTM D2622 & D4294	
ASI-SRO10-06-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	2,5	%	ASTM D2622 & D4294	
ASI-SRO10-07-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	3	%	ASTM D2622 & D4294	
ASI-SRO10-08-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	3,5	%	ASTM D2622 & D4294	
ASI-SRO10-09-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	4	%	ASTM D2622 & D4294	
ASI-SRO10-10-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	5	%	ASTM D2622 & D4294	