

Sulfur (S) (as Di-N-butyl sulfide) in Isooctane. Conc. Range: 0 - 0.300 wt%, 13 standards per set, 100 mL each (Standard for XRF)

Art. ID ASI-S(BS)I13-Set(13x100ML)
Unit 13 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-S(BS)I13-01-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0	%	ASTM D6334 & D2622	
ASI-S(BS)I13-02-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,001	%	ASTM D6334 & D2622	
ASI-S(BS)I13-03-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,002	%	ASTM D6334 & D2622	
ASI-S(BS)I13-04-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,003	%	ASTM D6334 & D2622	
ASI-S(BS)I13-05-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,005	%	ASTM D6334 & D2622	
ASI-S(BS)I13-06-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,01	%	ASTM D6334 & D2622	
ASI-S(BS)I13-07-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,02	%	ASTM D6334 & D2622	
ASI-S(BS)I13-08-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,03	%	ASTM D6334 & D2622	
ASI-S(BS)I13-09-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,04	%	ASTM D6334 & D2622	
ASI-S(BS)I13-10-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,06	%	ASTM D6334 & D2622	
ASI-S(BS)I13-11-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,1	%	ASTM D6334 & D2622	
ASI-S(BS)I13-12-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,2	%	ASTM D6334 & D2622	
ASI-S(BS)I13-13-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,3	%	ASTM D6334 & D2622	