

**Sulfur (S) (as Di-N-butyl sulfide) in Isooctane - Toluene mixture. Conc. Range: 0 - 0.1000 wt%, 10 standards per set, 100 mL each (Standard for XRF)**

Art. ID ASI-S(BS)IT10-Set(10x100ML)  
Unit 10 x 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-S(BS)IT10-01-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0	%	ASTM D6334 & D6445	
ASI-S(BS)IT10-02-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,0005	%	ASTM D6334 & D6445	
ASI-S(BS)IT10-03-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,001	%	ASTM D6334 & D6445	
ASI-S(BS)IT10-04-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,003	%	ASTM D6334 & D6445	
ASI-S(BS)IT10-05-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,005	%	ASTM D6334 & D6445	
ASI-S(BS)IT10-06-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,0075	%	ASTM D6334 & D6445	
ASI-S(BS)IT10-07-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,01	%	ASTM D6334 & D6445	
ASI-S(BS)IT10-08-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,03	%	ASTM D6334 & D6445	
ASI-S(BS)IT10-09-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,05	%	ASTM D6334 & D6445	
ASI-S(BS)IT10-10-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,1	%	ASTM D6334 & D6445	