

Sulfur (S) (as Di-N-butyl sulfide) in Isooctane. Conc. Range: 0 - 1000.00 ng/μL, 11 standards per set, 100ml each (Standard for XRF)

Art. ID ASI-SIO11-Set(11x100ML)
Unit 11 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-SIO11-01-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0	ng/uL	ASTM D5453	
ASI-SIO11-02-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,5	ng/uL	ASTM D5453	
ASI-SIO11-03-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	1	ng/uL	ASTM D5453	
ASI-SIO11-04-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	2,5	ng/uL	ASTM D5453	
ASI-SIO11-05-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	5	ng/uL	ASTM D5453	
ASI-SIO11-06-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	10	ng/uL	ASTM D5453	
ASI-SIO11-07-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	25	ng/uL	ASTM D5453	
ASI-SIO11-08-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	50	ng/uL	ASTM D5453	
ASI-SIO11-09-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	100	ng/uL	ASTM D5453	
ASI-SIO11-10-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	500	ng/uL	ASTM D5453	
ASI-SIO11-11-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	1000	ng/uL	ASTM D5453	