

**Sulfur (S) (as Di-N-butyl sulfide) in Isooctane. Conc. Range: 100.00 - 1000.00 ng/uL, 3 standards per set, 100 mL each (Standard for XRF)**

Art. ID ASI-SIO3(H)-Set(3x100ML)  
Unit 3 x 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-SIO3(H)-01-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	100	ng/uL	ASTM D5453	
ASI-SIO3(H)-02-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	500	ng/uL	ASTM D5453	
ASI-SIO3(H)-03-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	1000	ng/uL	ASTM D5453	