

Sulfur (S) (as Di-N-butyl sulfide) in Isooctane. Conc. Range: 5.00 - 100.00 ng/uL, 4 standards per set, 100ml each (Standard for XRF)

Art. ID ASI-SIO4(M)-Set(4x100ML)
Unit 4 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-SIO4(M)-01-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	5	ng/uL	ASTM D5453	
ASI-SIO4(M)-02-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	25	ng/uL	ASTM D5453	
ASI-SIO4(M)-03-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	50	ng/uL	ASTM D5453	
ASI-SIO4(M)-04-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	100	ng/uL	ASTM D5453	