

**Sulfur (S) (as Di-N-butyl sulfide) in Isooctane. Conc. Range: 0.0 - 10.00 ng/uL, 6 standards per set, 100ml each (Standard for XRF)**

Art. ID ASI-SIO6(L)-Set(6x100ML)  
Unit 6 x 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-SIO6(L)-01-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0	ng/uL	ASTM D5453	
ASI-SIO6(L)-02-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0,5	ng/uL	ASTM D5453	
ASI-SIO6(L)-03-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	1	ng/uL	ASTM D5453	
ASI-SIO6(L)-04-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	2,5	ng/uL	ASTM D5453	
ASI-SIO6(L)-05-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	5	ng/uL	ASTM D5453	
ASI-SIO6(L)-06-100ML	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	10	ng/uL	ASTM D5453	