

**Sulfur in Mineral Oil, Set of 6 calibration standards, concentrations of 0, 10, 50, 100, 200 and 500 ppm, 100 mL each (Standard for XRF)**

Art. ID ASI-SIN6C-Set(6x100ML)  
Unit 6 x 100 mL  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-SIN6C-01	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0	ppm	ASTM D7039	
ASI-SIN6C-02	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	10	ppm	ASTM D7039	
ASI-SIN6C-03	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	50	ppm	ASTM D7039	
ASI-SIN6C-04	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	100	ppm	ASTM D7039	
ASI-SIN6C-05	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	200	ppm	ASTM D7039	
ASI-SIN6C-06	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	500	ppm	ASTM D7039	