

Sulfur with Ultra Low Diesel Fuel, Set of 8 calibration standards, concentrations of 0, 5, 10, 15, 20 and 25ppm, 100ml each (Standard for XRF)

Art. ID ASI-ULSD0-25-Set(6x100ML)
Unit 6 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1202 III E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-ULSD0-25-01	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	0	ppm	EPA 80.580 to 80.585 and ASTM D7212, D7220	
ASI-ULSD0-25-02	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	5	ppm	EPA 80.580 to 80.585 and ASTM D7212, D7220	
ASI-ULSD0-25-03	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	10	ppm	EPA 80.580 to 80.585 and ASTM D7212, D7220	
ASI-ULSD0-25-04	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	15	ppm	EPA 80.580 to 80.585 and ASTM D7212, D7220	
ASI-ULSD0-25-05	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	20	ppm	EPA 80.580 to 80.585 and ASTM D7212, D7220	
ASI-ULSD0-25-06	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	25	ppm	EPA 80.580 to 80.585 and ASTM D7212, D7220	

